

DOWNLOADING IVEE AND VET-1 APPS




1

LOG INTO META ACCOUNT

- Log into your Meta Account on a desktop

OR

- Log into Meta Quest App on your phone



2

DOWNLOAD IVEE / VET-1

- IVEE link

 ivrytechnologies.com/ivee 

- VET-1 link

 ivrytechnologies.com/vet-1 


DOWNLOADING IVEE AND VET-1 APPS



3

ACKNOWLEDGE DISCLAIMER



- Accept the APP LABS disclaimer pop-up



4

CLICK "GET" BUTTON


- Click the blue "GET" button on the screen
- Allow time for download to process



DOWNLOADING IVEE AND VET-1 APPS

5

PUT HEADSET ON




6

CLICK THE APP LIBRARY



- Navigate to the app library on your headset
- Locate the IVEE and/or VET app
 - If you do not see the app it may still be finishing the download



DOWNLOADING IVEE AND VET-1 APPS

7

CLICK APP



- To install the app on your headset simply click the app
- Wait while the app installs. It may take up to 5 minutes

